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アクティニックLED630nm (290 $\mu\text{E}/\text{m}^2/\text{s}$)

Transient absorption changes induced by a continuous illumination with a 630 nm actinic LED (290 $\mu\text{E}/\text{m}^2/\text{s}$) following by their relaxation kinetics in the dark with a probing light at 520 nm.

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— P700

P700 oxidation induced by a continuous illumination with a 720 nm actinic LED (20 μ E/m²/s) followed by its reduction in the dark with a probing light peaking at 705 nm

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
Plastocyanin oxidation induced by a continuous illumination with a 720 nm actinic LED (20 μ E/m²/s) followed by its reduction in the dark with a probing light peaking at 740 nm

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A comparison of the plastocyanin and P700 time-course after normalization

Each curve is the result of a single sweep. The sample was preilluminated for 10 min.

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

Transient absorption changes induced by a continuous illumination with a 720 nm actinic LED (1230 μ E/m²/s) and a probing light at 705 nm (with a dark adapted spinach leaf)

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Transient absorption induced by a continuous illumination with a 720 nm actinic LED (120 μ E/m²/s) followed by its reduction in the dark : a detection at 820 nm for a dark adapted leaf (—) and for a light adapted leaf (preillumination of 3 minutes at 630 nm)(—)

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
Transient absorption changes induced by a continuous
illumination with a 720 nm actinic LED (190 μ E/m²/s)
and a probing light at 820 nm (—) and with a saturated
pulse at the end of the illumination (—)

The sample was preilluminated with a 630 nm actinic
light during 3 minutes

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
Fluorescence rise induced by a continuous illumination at
520 nm on a dark adapted leaf (Use of Fluo_59 board
accessory).

A few second of limiting intensity results in a large
increase of the fluorescence yield indicating the
progressive reduction of the plastoquinone pool.

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Fluorescence yield changes induced by a continuous
illumination at 520 nm and in the dark.

The maximum fluorescence yield is assessed by a
surimposition of short intense green light pulses (7900
 μ E/m²/s).

Access to Fo, Fm and Fm' values

Note: all Fm' values are not represented on the curve.

The associated Genty parameters Φ PSII* and the Non
Photochemical Quenching

Parameters NPQ* in function of the time can then be
calculated once Fo, Fm and Fm' are determined.

The opposite graph shows the evolution of ξ Φ PSII*
parameter in function of the time(—)


Φ PSII* was determined with the following equation:

$$\Phi PSII = \frac{Fm' - Fo}{Fm'}$$

NPQ evolution is as well followed in function of the
time (—)

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
Figure 1 : Transthylakoid pH variations

The measurement was done on young leaf spinach and shows a transient absorption changes induced by a continuous illumination with a 630 nm actinic LED, followed by the relaxation kinetic in the dark. The probing light wavelength is 520 nm. The relaxation signal in the dark (i.e. when the actinic light is switched-off) reflects the pH variation in the transthylakoid region.


Figure 2 : use of Cytochrome eukaryote kit

Measurements were done on young spinach leaves to follow the cyt b signal:

- Yellow curve (fair): Use of actinic orange light ($150 \mu E/m^2/s$) and detection at 546 nm
 - Green curve: Use of actinic orange light ($150 \mu E/m^2/s$) and detection at 554 nm
 - Light green curve: Use of actinic orange light ($150 \mu E/m^2/s$) and detection at 573 nm.
- The deconvolution of these 3 signals allows to obtain the pure kinetic of cytochrome f reflecting the redox changes (curve below)


Figure3: Cyt f: pure signal

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